



Since its founding, Ushio Medical Systems (formerly Nathaniel Group) has supplied complex mechanical and electronic equipment for the medical industry. Ushio Medical Systems has expanded its expertise to include optical and illumination systems. Theia™ series products combine an LED illuminator with an integrated video processing unit (VPU) for original equipment manufacturers (OEMs) in microendoscopy applications. The Theia 40 has an outside diameter of 1.1mm and has auto-LED brightness. The Theia 160 has auto-gain and auto-exposure features for endoscopes as small as 1.7mm outside diameter (OD). Both Theia models include a 40k-pixel resolution or 160k-pixel resolution system with the latest CMOS image sensor technology.

VPUs on each unit have a programmable measurement window or region of interest (ROI). The VPU overlays a virtual window border over an image allowing ROI changes to the sides of the center window border without impacting the image size or frame rate. (see Figure 1).

In general, the Theia VPU uses a defined target gray value and calculates the mean gray value within the ROI. The firmware then decides whether to increase the brightness, keep the brightness, or decrease the brightness with the optimal image setting.




Figure 1: Theia VPU with virtual window enabled. The ROI has no impact on image size or frame rate.




Figure 4: From left to right, Theia 40 Low, Default, and High auto LED brightness modes.

Theia 40 Auto-LED Brightness

The Theia 40 is ideal for the smallest endoscopes available, as small as 1.1mm OD. The small OD endoscope does not allow for auto-gain or auto-exposure, but the Theia VPU has an auto-LED brightness feature for ease of use.

Auto Mode

There are three different target values in auto mode: low, default and high. The user can toggle the Auto Mode and the Theia VPU will re-balance the LED illumination intensity to the selected target value (see Figure 4).

Manual Mode

In manual mode the user can select the corresponding LED illumination intensity based on the best image and the Theia 40 VPU will re-balance to the settings (see Figure 5).




Figure 5: From left to right, Theia 40 Low, Medium and High manual modes.

Theia 160 Auto-Gain & Auto-Exposure

Increase In brightness

The Theia 160 will first increase the exposure time. If the max exposure time is reached, then the gain will be increased by pixel voltage amplification.

Decrease In brightness

The Theia 160 will first decrease the gain that lowers image noise. If the minimum gain is reached, then the exposure time is lowered.

Auto Mode

There are three different target values in auto mode: Low, Default and High. The user can toggle the Auto Mode and the Theia VPU will re-balance the mean to the selected target gray value (see Figure 2).




Figure 2: From left to right, Theia 160 Low, Default and High auto modes.

Manual Mode

There are also three fixed exposure and gain settings in manual mode. Users can select the corresponding exposure and gain value while the Theia VPU will re-balance to the settings. Compare Auto and Manual modes to see how settings give the user more flexibility for their application (see Figure 3).




Figure 3: From left to right, Theia 160 Low, Medium, and High Manual modes.

Summary

Both the Theia 40 and Theia 160 offer the highest video quality and resolution with the smallest OD endoscope sizes available on the market. The standard Theia VPU features auto-gain and auto-exposure, or auto-LED brightness depending on the model required for the end-user's application.

Custom VPU, LED, and videoscope options are also available, including trade dress for OEMs. Let us know how Ushio can help you find the right solution.

Ushio America, Inc.
Author: Brandon Morioka

For Further Information - Please Contact
Keith Harrison - Assistant General Manager
T: +1.802.877.2182
kharrison@ushio.com
www.ushio.com/medical-systems